
ALBORA Red Light Belt for Dogs

User Manual
Model VL-MAT-S01



Limited Warranty
Page 12

BEFORE YOU START

Thank you for choosing ALBORA. Please read this manual before first use and keep it for future reference.

In the box

- ALBORA LED belt
- 800 mm (31.5 in) strap
- 1000 mm (39.4 in) strap
- USB-C charging cable (1 m)
- Storage pouch
- This manual
- Session & Wellness Tracker
- Thank You & Support Card


No wall adapter is included.
See page 6 for a compatible adapter.

The two straps are alternatives. Use one strap at a time — choose the one that fits your dog (page 4).

Quick start

1. Charge the belt fully before first use (page 6).
2. Choose and thread the strap that fits your dog (page 4).
3. Wrap the belt around your dog with the LED side against the body (page 5), or lay it flat as a mat (page 8).
4. Press the button to start. The belt switches off automatically after 15 minutes. Stay with your dog throughout.

KNOW YOUR BELT




1. LED side — Place this side against your dog's body.
2. Battery end — This end of the belt houses the battery and has no LEDs.
3. Slots — Two reinforced slots for attaching the strap. Use one slot at a time.
4. Power button — Turns the belt on, changes modes and shows the battery status.
5. USB-C port — For charging.
6. Strap — Has a smooth side and a fuzzy side. The grey X marks the end you thread first.

Belt size: 46 × 21.6 cm (18.1 × 8.5 in)

CHOOSE AND THREAD THE STRAP

1. Measure your dog




Measure around the **widest part of the chest, just behind the front legs.**

Strap	Approx. chest girth
800 mm	16.5 – 24.5 in (42 – 62 cm)
1000 mm	23.5 – 31.5 in (60 – 80 cm)

Fitting ranges are approximate and may vary depending on your dog's body shape and coat. Choose by measurement, not by breed or weight. If your dog is between sizes, try the 800 mm strap first.


2. Thread the strap



1. **Start with the grey X.** Find the end of the strap marked with the grey X. It is the only grey mark on the strap.
2. **With the outer side facing up, feed this end through the slot next to the power button.**
3. **Fold the end back onto the strap.** Press the grey X firmly onto the fuzzy side to form a small loop around the slot. **The red X faces out.**
4. **Gently pull the strap** to check that the fastener holds securely.


Leave the other end free.
You will close it once the belt is on your dog (page 5).
Watch the setup video: getalbora.com/setup

PUT IT ON YOUR DOG



1. Place the **LED side flat against your dog's body** over the area you want to cover.
2. Bring the strap around **over the outside of the belt** — never between the LEDs and your dog's skin.
3. **Close the strap:** press the free end onto the fuzzy side of the strap. Adjust it so it is **secure and comfortable — neither loose nor tight.** Once fastened, only red X marks show.
4. **Two-finger check:** you should be able to slide two fingers under the strap. It must not restrict breathing or movement.

Smaller dogs: if the ends of the belt overlap, place the battery end (the end without LEDs, see page 3) on the outside. Keep the overlap within the section without LEDs.



If thick fur prevents the belt from lying flat, ask your veterinarian for advice. Do not shave your dog or press the belt down to improve contact.

CHARGING

1. Make sure the belt is switched off and the USB-C port and cable are **clean and dry.**
2. Connect the **included USB-C cable** (recommended) to the belt and to a **USB power adapter or power bank with 5 V output (2 A recommended), certified by UL or ETL.**
3. While charging, the button light **flashes red**. When it turns **steady red,** the belt is fully charged.
4. A full charge takes about **8 hours.** Charging stops automatically when the battery is full. Disconnect the cable once the light is steady red.

Before first use: charge the belt fully.

Battery low: if the button light flashes red while you are using the belt (not connected to a charger), the battery is running low. Recharge it.

Belt input: 5 V⎓ 1 A, USB-C.

⚠

WARNING — Lithium-ion battery

- Never charge while your dog is wearing the belt.
- Do not charge near water, under bedding, or on soft surfaces that trap heat.
- Charge on a hard, stable surface. Do not leave the belt charging unattended or while you sleep.
- Do not use or charge the belt if it is swollen, leaking, unusually hot or damaged.
- Do not open the belt or try to replace the battery.

CONTROLS AND LIGHT MODES

The belt has one button and two light modes.

Steady light — short press from off. Button light: **green.**
Pulsed light — with the belt switched on in any light setting, press and hold the button for **about 3 seconds.** Button light: **blue.**

To return to steady light: switch the belt off using short presses, then short-press the button to switch it on again.

Switching on starts with red + near-infrared light. While the belt is on in either mode, each short press advances to the next step in this sequence:

Step	Light
1	Red + near-infrared (660 nm + 850 nm)
2	Red only (660 nm)
3	Near-infrared only (850 nm)
4	Off

To stop a session early: use short presses to advance through the remaining settings until the belt switches off. The number of presses depends on the current setting. Holding the button activates pulsed light; it does not switch the belt off.

The belt switches off automatically **15 minutes after you switch it on.** Changing the light setting does not restart the timer.

Button light

Button light	Meaning
Green	Steady-light mode
Blue	Pulsed-light mode
Flashing red, connected to a charger	Charging
Steady red, connected to a charger	Fully charged
Flashing red during use	Low battery — recharge


Near-infrared light is invisible to the human eye. In near-infrared-only mode, the LEDs may look off, **but the button indicator remains illuminated while the belt is operating. Do not look into the LEDs to check.**

SESSIONS AND SUPERVISION

- For dogs only. An adult must stay with the dog for the whole session.
- Let your dog get used to wearing the belt switched off before the first session.
- Stop the session if your dog tries to remove the belt or shows any sign of discomfort.
- The belt switches off automatically after 15 minutes. You can stop it earlier if needed (page 7). Do not restart the belt to extend a session unless your veterinarian advises it.
- The automatic 15-minute shut-off does not determine the right session length for your dog. Follow your veterinarian's advice on session duration and frequency.
- Ask your veterinarian whether the belt is suitable for your dog and how often to use it, especially if your dog has a health condition or is receiving treatment.

Use it as a mat

You can also lay the belt flat and let your dog lie on it.



1. Remove the strap.
2. Place the belt on a flat surface, such as the floor, a sofa or a bed, with the **LED side up.**
3. Switch it on and let your dog lie with the area you want to cover resting on the LEDs.
4. Keep the **battery end** (the end without LEDs) clear of your dog's body.

- For dogs up to **110 lb (50 kg),** lying or resting only. Do not let your dog walk, stand, jump or play on the belt. Do not stand, kneel or sit on it yourself.
- Do not cover the belt or your dog with a blanket during the session.
- Do not use it as a mat while it is charging.
- Stop the session if your dog chews, scratches or digs at the belt.

ALBORA does not replace veterinary care.

SAFETY

⚠ **WARNING** — Eyes. Never shine the LEDs into human or animal eyes. Never stare into the LEDs.

⚠ **WARNING** — Battery. Do not puncture, crush, heat or burn the belt. Do not expose it to fire or high temperatures.

⚠ **CAUTION** — Health. Ask your veterinarian before use if your dog is pregnant, has a suspected tumor, takes medication that increases light sensitivity, or has any other medical condition.

⚠ **CAUTION** — Skin. Do not use over broken skin, open wounds or a recent surgical site unless your veterinarian has advised it. Stop use if you notice redness, heat, irritation or unusual behavior.

⚠ **CAUTION** — Supervision. Never leave the belt on an unattended dog. Keep the straps and cable away from chewing and tangling.

NOTICE — Product care. Keep the belt dry. Do not use it if it is damaged. For dogs only. Not a toy. Keep out of reach of children.

CARE, STORAGE AND TROUBLESHOOTING

Cleaning

1. Switch the belt off and **disconnect the charging cable.**
2. **Fabric parts and straps:** clean gently with a damp towel or a soft toothbrush.
3. **Outer surface and LED side:** gently wipe them with a wet wipe only. Do not clean the LEDs directly or try to remove the protective surface over them.
4. **Do not submerge the belt** in water or put it in a washing machine. Do not spray the belt. Keep moisture out of the power button and USB-C port. Do not use cleaning products or solvents.
5. Let the belt **air-dry** out of direct sunlight. Make sure it is completely dry before using or charging.

Storage



Roll — do not fold. Folding the belt can disconnect the LEDs from the circuit. The belt can curve gently around your dog, but never crease it. Store it rolled in a cool, dry place, out of reach of pets and children.

Troubleshooting

Problem	What to do
Belt won't turn on	Charge it fully. Check the cable and power adapter.
Belt switches off during a session	If about 15 minutes have passed, the automatic timer has ended the session. If it stops earlier, charge the belt. If the problem continues, contact support.
Button light doesn't come on when charging	Check that the port and cable are clean, dry and undamaged. Try another 5 V adapter.
LEDs look off in near-infrared mode	Normal — near-infrared light is invisible. In near-infrared-only mode, the button indicator remains illuminated even though the light is not visible.
Strap loosens or rubs	Stop the session. Refit it, or switch to the other strap.

Still not working? Stop using the belt and contact us at support@getalbora.com.

SPECIFICATIONS AND FCC INFORMATION

Model	VL-MAT-S01
Belt size	46 × 21.6 cm (18.1 × 8.5 in)
LEDs	120 LEDs, each combining one red chip (660 nm) and two near-infrared chips (850 nm)
Irradiance	70 mW/cm² at the LED surface, through the protective cover (red + near-infrared combined)
Battery	Built-in rechargeable lithium-ion, 3.7 V, 3600 mAh, 13.32 Wh. Not user-replaceable
Charging input	5 V⎓ 1 A, USB-C
Auto-off	15 minutes
Straps	800 mm and 1000 mm
Weight	9.8 oz (278 g), belt without straps
Made in	China

FCC Supplier's Declaration of Conformity

47 CFR §2.1077 Compliance Information
Unique identifier:
ALBORA Red Light Belt for Dogs, Model VL-MAT-S01

Responsible party (U.S.):
Ventaris Store LLC
412 W 7th St, Clovis, NM 88101, United States

Contact:
support@getalbora.com

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

NOTE: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:
- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

WARRANTY, SUPPORT AND RECYCLING

24-Month Limited Warranty

Warrantor. This warranty is provided by Ventaris Store LLC, 412 W 7th St, Clovis, NM 88101, United States · support@getalbora.com.

Who is covered. This warranty covers the original purchaser of a new ALBORA belt bought from Ventaris Store LLC or an authorized seller in the United States.

What is covered. Defects in materials and workmanship under normal use for **24 months from the date of purchase.** A replacement belt is covered for the remainder of the original warranty period.

What is not covered. Normal wear of the straps and hook-and-loop fasteners, or damage caused by misuse, accidents, chewing or scratching, sharp folds or creases contrary to this manual, liquid ingress or immersion, incompatible charging equipment, or improper repairs or modifications. Cleaning the protective surfaces as instructed in this manual is permitted. Using a compatible third-party cable or obtaining third-party service does not, by itself, void this warranty; damage caused by such equipment or service is not covered.

What we will do. We will replace the belt or refund the purchase price, at our option. For claims covered by this warranty there is **no charge to you:** we pay the shipping of any replacement within the United States. If we need the defective belt back, we will send you a prepaid shipping label. **Please do not return the belt until you have received instructions from us.**

How to make a claim. Email **support@getalbora.com** with your order number, a short description of the problem and, where possible, a photo or short video. We will send you an initial reply within 2 business days with next steps. You will need proof of purchase.

Implied warranties. Any implied warranties, including the implied warranties of merchantability and fitness for a particular purpose, are limited to the duration of this written warranty. Some states do not allow limitations on how long an implied warranty lasts, so the above limitation may not apply to you.

Incidental and consequential damages. To the extent permitted by applicable law, we are not liable for incidental or consequential damages. Some states do not allow the exclusion or limitation of incidental or consequential damages, so the above limitation or exclusion may not apply to you.

This warranty gives you specific legal rights, and you may also have other rights which vary from state to state.

Full warranty terms: getalbora.com/warranty
To request a free printed copy, write to Ventaris Store LLC, 412 W 7th St, Clovis, NM 88101, United States.

Recycling

This product contains a lithium-ion battery. Do not dispose of it in household waste. Take it to a battery or electronics recycling point. Find a drop-off location at call2recycle.org.

ALBORA is a brand of Ventaris Store LLC · 412 W 7th St, Clovis, NM 88101, United States · support@getalbora.com · Made in China